

Revive

HexCell rPET Trays



CODE	DESCRIPTION	COLOUR	DIMENSIONS	UOM/QTY
REvhxc75	Revive HexCell rPET Tray	Clear	7 x 5"	CTN/680
REvhxc85	Revive HexCell rPET Tray	Clear	8 x 5"	CTN/660
REvhxc87	Revive HexCell rPET Tray	Clear	8 x 7"	CTN/500
REvhxc115	Revive HexCell rPET Tray	Clear	11 x 5"	CTN/420
REvhxc119	Revive HexCell rPET Tray	Clear	11 x 9"	CTN/280
REvhxc1310	Revive HexCell rPET Tray	Clear	13 x 10"	CTN/260
REvhxc1411	Revive HexCell rPET Tray	Clear	14 x 11"	CTN/170

*All sizes now available with rolled edge feature.



Innovative Liquid Retention Technology

Revive HexCell Trays are designed to capture and retain liquid in the cells at the base of the tray - helping to prevent leakage through the overwrap and eliminating the need for soaker pads.



Recycled Content

Contains Recycled Content

Revive HexCell Trays are made using recycled materials rather than using only virgin material. Purchasing products made from recycled materials helps to support the recycling economy.



Widely Recyclable

Once rinsed thoroughly, Revive HexCell Trays can be disposed of in co-mingled recycling. This means they can be collected, sorted and recycled successfully in practice and at scale in Australia and New Zealand.



Exclusively Distributed by Bunzl

www.bunzl.com.au/revive

